
# THE LOCKER™

Screw Augmentation System

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# THE LOCKER™ Screw Augmentation System

EN

## INSTRUCTIONS FOR USE *For the attention of the surgeon*

### **PRESENTATION**

The Locker is a sterile disposable device for the preparation and application of acrylic resin in screw augmentation. The device is contained in a hermetically sealed blister, sterilised with ethylene oxide. The phial liquid is sterilised by means of filtration; the needle is enclosed in a sealed double container, sterilised with ethylene oxide. The system enables the injection of acrylic resin into the bone through a screw. The injection needle has a 13 gauge diameter and has to be used with an appropriate screw.

The Locker consists of:

- Radio-opaque, high-fluidity acrylic resin
  - \* 20 g container of sterile powder
    - Polymethylmethacrylate 67,5% w/w
    - Barium sulphate 30,0% w/w
    - Benzoyl peroxide 2,5% w/w
  - \* 9,4 g phial of sterile liquid
    - Methylmethacrylate 99,1% w/w
    - N,N-dimethyl-p-toluidine 0,9% w/w
    - Hydroquinone 75 ppm
- Injection system:
  - Screw Augmentation gun
  - Aspiration tube
  - Connection tube
  - Screw augmentation needle

The injection system uses Luer-lok type connections.

### **INTENDED USE**

The Locker is a system intended to increase the fixation strength of screws to the bone (screw augmentation).

### **CONTRAINDICATIONS**

Active infections  
Emorrhagic diathesis  
Fractured or partially resected vertebral bodies  
Fever

The use of The Locker must be carefully considered in presence of hypersensitivity to monomer or to any of the other components of the resin.

### **USE DURING PREGNANCY, BREASTFEEDING AND IN CHILDREN**

There are no tests which demonstrate the utilisation safety of The Locker during pregnancy or breast-feeding. The Locker should not be used in the first three months of pregnancy; for the remaining gestation period, it should only be used in life-endangering situations.

The use in children is advised only when it is believed impossible to operate through other forms of intervention.

### **PRECAUTIONS FOR USE**

- Read this instruction booklet very carefully.
- Read the instructions for use of the appropriate implant system.
- Make sure that the inner package is undamaged and that the components are integral, i.e. the powder should not present yellow or brown discolouring and the liquid should not be syrupy.
- Temperature has a major effect on the preparation characteristics of any acrylic resin. Temperatures of more than 23°C for the product, the preparation accessories and the environment accelerate the various stages in the preparation procedure. Lower temperatures retard the preparation stages. Before using The Locker it is strongly advised to make sure that the package is stored at a temperature of 23°C ± 1°C for the previous 24 hours.
- Excessive exposure to a high humidity environment (>70%) may lead to an increase in viscosity, thus speeding up resin preparation and application times.
- Do not open the phial of liquid over the mixing bowl to avoid the risk of glass fragments entering the dough.
- NEVER add other substances or foreign bodies to the acrylic resin.

To ensure safe application of The Locker, surgeons must be aware of the properties of the product, its characteristics, applicational limitations and correct application technique. It is strongly recommended that the surgical team carry out practical trials prior to use in patients under the same instrumental and environmental conditions.

Clinical data demonstrate the need to maintain strictly aseptic surgical techniques. It is important to keep in mind that the possible deep infection of a surgical wound is a serious risk for the successful outcome of the surgical operation.

Prolonged respiration of the vapours may cause drowsiness; excessive exposure to monomer vapours may irritate the respiratory apparatus and the eyes. Avoid monomer contact with the skin and mucous membranes. Cases of contact dermatitis have been observed in susceptible subjects. It is therefore advisable to wear a second pair of surgical gloves and scrupulously observe the instructions for mixing the components in order to reduce the possibility of reactions caused by hypersensitivity. Do not allow the liquid component to come into direct contact with accessories in rubber or elastomers.

The vapours of the liquid component may affect soft contact lenses.

The liquid component is flammable and volatile and for this reason the operating theatre must be correctly ventilated. The liquid component and/or its vapours must never be directly exposed to naked flames or heated materials.

**Caution:** Acrylic resins reach temperatures higher than physiological temperatures during the exothermic polymerisation reaction. Following surgery, if any form of infection should arise, patients must immediately consult their doctors to reduce the risk of infection.

### **TIMING FOR THE PREPARATION AND APPLICATION**

Times indicated in the table were obtained by performing tests under temperature controlled conditions of 23°C ± 1°C and a relative humidity 50% ± 10%.

Any increase in temperature of the working environment, of the resin components and of the preparation accessories over 23°C reduces the times reported in the table. Equally, lower temperatures increase such times.

<b>Operations</b>	Mixing	Transfer into the Gun	Waiting	Application
<b>Duration</b>	60''	60''	60''÷180''	600''

**Warning:** Start the application when the resin has a pasty consistency in order to avoid cement leakage due to low viscosity.

During the application of the resin, the viscosity of the resin can be fitted to the implant site by anticipating or delaying the moment of application. Do not exceed the time period indicated, since the fluidity of the resin decreases over time, and it becomes more difficult to dispense.

#### **INSTRUCTIONS FOR USE**

• **Resin preparation:** Set the container of powder on a flat surface and open the cap. Open the vial and pour all of the liquid into the powder. Close the container and manually shake it horizontally until a liquid and homogenous mass is achieved. By overturning the container vertically, check if any of the mass is packed at the bottom of the container. In this case, shake energetically for obtaining an homogeneous mass. Set the container down again, remove the cap.

**Caution:** Never arbitrarily modify the ratios between the liquid and solid components.

• **Loading the Screw Augmentation Gun:** grip the Gun, attach the aspiration tube and immerge the tip of the tube into the container with the resin in a fluid state. Squeeze the Gun handle firmly (thus releasing the safety catch) and pull the piston upwards until the Gun is completely full. Release pressure on the handle.

**Caution:** A subsequent reloading of the Gun is not possible.

• **Bleeding the system:** For bleeding the filled system the tip of the Gun has to point upward. Turn the threaded piston until cement exits the tip.

• **Removal of the aspiration tube:** Turn the ring nut on the tip of the Gun to free the Luer-lok, and then pull out the aspiration tube.

• **Optional:** attachment of the connection tube connecting the Gun and the injection needle through the Luer-lok connection.

The correct choice of surgical technique will depend on the type of screw used:

1) **Use with non-cannulated synthesis devices:** remove the stylet from the needle and connect the connection tube to the needle, tightening the ring nut on the tube to the needle. Inject the resin directly into the bone cavity before positioning the screw (see example 1).

2) **Use with cannulated synthesis devices equipped with Luer-lok connections:** attach the connection tube to the screw head through the Luer connection and inject the resin. Check that the connection between the screw and the connection tube is firmly in place before beginning the injection (see example 2).

3) **Use with cannulated synthesis devices without Luer-lok connections:** place the injection needle with the stylet directly into the screw and carefully check that the cannula needle does not contain any foreign bodies. Remove the stylet from the needle and connect the connection tube to the needle, tightening the ring nut. Begin the resin injection (example 3).

**Warning:** before injecting the cement, ensure that the connection between the screw and the needle is firmly fixed.

• **Supplying resin:** Resin must be supplied to the bone by turning the threaded piston clockwise and continuously.

**Warning:** Monitor the filling procedure using high-quality fluoroscopy. The injected cement volume is related to the screw used. Carefully read the instructions for use of the appropriate implant system.

• **Removal of the injection needle** before resin polymerisation

#### **WARNINGS**

• Maintain at a temperature below 25°C, with a relative humidity below 70%, and away from the light.

• Sterility is assured only if the unit container is not damaged or opened.

• Do not re-sterilise the set or any of the components.

• Do not use the product if the powder has a yellowish or brownish colour or if the liquid is syrupy. These two conditions indicate that the product has not been stored correctly.

• The use of different injection systems must be carefully evaluated by the surgeon.

• The screw augmentation procedure must be carefully evaluated because the resin can only be removed with a surgical operation.

#### **ADVERSE EVENTS**

The following negative reactions may arise with the use of the acrylic resin, but they are not directly due to the acrylic resin itself.

The surgeon must be aware of these possible events and be ready to treat them if they appear.

• cardiac arrest

• myocardial infarct

• pulmonary embolism

• temporary drop in blood pressure

• haemorrhage and haematoma

• infection of surface/deep surgical lesion

• short-term irregularity of the cardiac rhythm

• sudden death

• diffusion of the resin outside the vertebral body: in the perivertebral veins (pulmonary embolism), in the epidural plex (mielopathy, radicolopathy), in the intervertebral disc.

• ischemic necrosis of the femoral head

#### **WASTE DISPOSAL**

Disposal of the device and its components must be carried out in compliance with local waste-disposal regulations.

Before elimination, allow the excess resin to polymerize.

# THE LOCKER™ Screw Augmentation System

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TECRES S.p.A.