

Series 7000 Large Power System from MicroAire®



PNEUMATIC LARGE POWER STERILIZATION CASE 9700-000

Features Durable polymer and aluminum construction

Base holds 3 handpieces and hose

Insert tray holds 11 couplers and pin driver

10" x 20" x 4" (25cm x 50cm x 10cm)

AIR HOSES & ACCESSORIES

MicroAire offers a wide variety of Air Hoses and accessories for use with your Series 7000 instrument, and other pneumatic instrument systems.

For more information, please see your MicroAire Catalog, contact your MicroAire Sales Representative, or visit us online at www.microaire.com

Instruments & Accessories



Catalog Number	Description
7105-200	Drill/Reamer with Hall® Connector
7205-200	Oscillating Saw with Hall® Connector
7405-200	Reciprocating Saw with Hall® Connector

Handpieces are also available with the following air hose connectors:
Linvatec/3M® "Maxi", A.O. Synthes®, Stryker®

7013-000	"Whisper" Hall® Style Air Hose
9013-000	Hall® Style Air Hose
7100-001	Zimmer® Reamer Coupler (Hudson® Modified Trinkle), 250 RPM
7100-002	Trinkle Coupler, 750 RPM
7100-003	Jacobs Coupler 1/4", 750 RPM
7100-004	Keyless Coupler 1/4", 750 RPM
7100-005	Keyless Coupler 1/8", 750 RPM
7100-006	Synthes® Quick Connect Coupler, 750 RPM
7100-008	Hudson® Reamer Coupler, 250 RPM
7100-010	Jacobs Coupler 1/4", 250 RPM
7100-011	Harris/Aesculap® Reamer, 250 RPM
7100-011	Synthes® Reamer Coupler, 250 RPM
7100-112	Synthes® Non-Reversing Reamer Coupler, 250 RPM
7100-045	Wire/Pin Driver Coupler, 1.6mm- 2.8mm
7100-050	Pin Driver Coupler, 3.2mm- 4.0mm
7100-051	Zimmer® Coupler (Hudson® Modified Trinkle), 750 RPM
7100-058	Hudson® Coupler, 750 RPM
7100-474	Pulse Lavage Coupler
9700-000	Sterilization Case



MicroAire Surgical Instruments LLC
1641 Edlich Drive
Charlottesville, Virginia 22911
U.S.A.

(800) 722 0822
(434) 975 8000
Fax: (800) 648 4309
Fax: (434) 975 4131
www.microaire.com



Product Information

SERIES7000™
Pneumatic Power





7105-200 Drill/Reamer

DRILL

Speed **750 RPM**
Torque **50 in.lbs (5.6 Nm)**

REAMER

Speed **250 RPM**
Torque **150 in.lbs (16.9 Nm)**

Drilling



Keyless Coupler 1/4"
Speed **750 RPM**
Catalog # **7100-004**



Keyless Coupler 1/8"
Speed **750 RPM**
Catalog # **7100-005**



Synthes® Quick Connect Coupler
Speed **750 RPM**
Catalog # **7100-006**



Zimmer® (Hudson®) Modified Trinkle Coupler
Speed **750 RPM**
Catalog # **7100-051**



Hudson® Coupler
Speed **750 RPM**
Catalog # **7100-058**



Trinkle Coupler
Speed **750 RPM**
Catalog # **7100-002**



7205-200 Oscillating Saw

SAW

Speed **13,000 CPM**
Arc **4.5°**
12 Position Rotating Head

Reaming



Zimmer® (Hudson®) Modified Trinkle Reamer Coupler
Speed **250 RPM**
Catalog # **7100-001**



Hudson® Reamer Coupler
Speed **250 RPM**
Catalog # **7100-008**



Jacobs® Coupler 1/4"
Speed **250 RPM**
Catalog # **7100-010**



Harris/Aesculap® Reamer Coupler
Speed **250 RPM**
Catalog # **7100-011**



Synthes® Reamer Coupler
Speed **250 RPM**
Catalog # **7100-111**



Synthes® Non-Reversing Reamer Coupler
Speed **250 RPM**
Catalog # **7100-112**



7405-200 Reciprocating Saw

SAW

Speed **14,000 CPM**
Blade Stroke **3.2mm**
360° Blade Rotation

Pulse Lavage



Pulse Lavage Coupler
Catalog # **7100-474**
Uses 5740 tubing sets

Pin Driving



Jacobs® Coupler 1/4"
Speed **750 RPM**
Catalog # **7100-003**



Wire/Pin Driver Coupler
Speed **750 RPM**
Catalog # **7100-045**



Pin Driver Coupler
Speed **750 RPM**
Catalog # **7100-050**



Precision, Quality, MicroAire