

MICROAIRE

Introducing the new large power line from MicroAire. Our new handpieces function as both cordless AND electric instruments, offering the ultimate in flexibility. It's the versatility the industry has always wanted with the security of the MicroAire name.

Direct Power Console 7025

- Features:**
- Simple dedicated electric power supply
 - Autoclavable Console Pack
 - Backup for battery



Universal Battery Charger 7500-700

- Features:**
- 14.4 Volt Nickel Metal Hydride (NiMH) rechargeable, autoclavable battery
 - UL 2601-1 Compliant battery charger with four charging bays
 - NiMH batteries have more power and longer run time than conventional NiCd batteries
 - Also charges 6640-710 SmartDriver batteries



INSTRUMENTS & ACCESSORIES

CATALOG NUMBER	DESCRIPTION
7025	Direct Power Console (DPC)
7007-000	DPC Console Pack
7100-001	Zimmer Reamer Coupler, 250 RPM
7100-002	Trinkle Coupler, 750 RPM
7100-003	Jacobs Coupler, 750 RPM
7100-004	Keyless Coupler 1/4", 750 RPM
7100-005	Keyless Coupler 1/8", 750 RPM
7100-006	A.O. Synthes Coupler, 750 RPM
7100-008	Hudson Reamer Coupler, 250 RPM
7100-010	Jacobs Coupler, 250 RPM
7100-011	Harris/Aesculap Reamer, 250 RPM
7100-111	A.O. Synthes Reamer Coupler, 250 RPM
7100-112	A.O. Synthes Non-Reversing Reamer Coupler, 250 RPM
7100-045	Pin Driver Coupler, 1.6mm- 2.8mm
7100-050	Pin Driver Coupler, 3.2mm- 4.0mm
7100-051	Zimmer Coupler, 750 RPM
7100-058	Hudson Coupler, 750 RPM
7100-474	Pulse Lavage Coupler
7500	Drill/Reamer
7500-700	Universal Battery Charger
7500-710	Large Power Battery Pack
7600	Oscillating Saw
7700	Reciprocating Saw
4100-021	Automatic Screwdriver Bit Holder for 4100-024/025/026
4100-023	Automatic Screwdriver Bit Holder for 4100-027
4100-024	Screwdriver Bit, Single Slot
4100-025	Screwdriver Bit, Phillips
4100-026	Screwdriver Bit, Cruciate
4100-027	Screwdriver Bit A.S.I.F. Trinkle Shank
9720-000	Sterilization Case

1641 Edlich Drive Charlottesville, VA 22911
 TOLL FREE (800) 722-0822 (434) 975-8000
 FAX (800) 648-4309 (434) 975-4131
 Visit us online at: www.microaire.com



INTRODUCING
SERIES 7000

Cordless. Electric. We offer both in one instrument.





7500 DRILL/REAMER

DRILL
 SPEED 750 RPM
 TORQUE 50 in. lbs (5.6 Nm)

REAMER
 SPEED 250 RPM
 TORQUE 150 in. lbs (16.9 Nm)



7600 OSCILLATING SAW

SAW
 SPEED 15,000 CPM
 ARC 4°

12 POSITION
 ROTATING HEAD



7700 RECIPROCATING SAW

SAW
 SPEED 15,000 CPM
 BLADE STROKE 3.2mm

360° Blade Rotation

DRILLING



Jacobs Coupler 750 RPM 7100-003
NAME SPEED CATALOG #



Zimmer Coupler 750 RPM 7100-051
NAME SPEED CATALOG #



Keyless Coupler 1/4" 750 RPM 7100-004
NAME SPEED CATALOG #



Hudson Coupler 750 RPM 7100-058
NAME SPEED CATALOG #



A.O. Synthes Coupler 750 RPM 7100-006
NAME SPEED CATALOG #



Keyless Coupler 1/8" 750 RPM 7100-005
NAME SPEED CATALOG #



Trinkle Coupler 750 RPM 7100-002
NAME SPEED CATALOG #

REAMING



Harris/Aesculap Reamer Coupler 250 RPM 7100-011
NAME SPEED CATALOG #



Synthes Non-Reversing Reamer Coupler 250 RPM 7100-112
NAME SPEED CATALOG #



Zimmer Reamer Coupler 250 RPM 7100-001
NAME SPEED CATALOG #



Hudson Reamer Coupler 250 RPM 7100-008
NAME SPEED CATALOG #



Synthes Reamer Coupler 250 RPM 7100-111
NAME SPEED CATALOG #



Jacobs Coupler 250 RPM 7100-010
NAME SPEED CATALOG #

PIN DRIVING



Jacobs Coupler 250 RPM 7100-003
NAME SPEED CATALOG #



Pin Driver Coupler 750 RPM 7100-045
NAME SPEED CATALOG #
 NOTE: ACCEPTS PINS FROM 1.6 MM TO 2.8MM



Pin Driver Coupler 750 RPM 7100-050
NAME SPEED CATALOG #
 NOTE: ACCEPTS PINS FROM 3.2MM TO 4.0MM

K-Wires & Steinmann Pins SIZE FROM 1.6MM TO 4.0MM

**SCREWDRIVING COUPLERS
 BITS & HOLDERS**



Trinkle Coupler 750 RPM 7100-002
NAME SPEED CATALOG #



Screwdriver Bit Holder 4100-021
NAME CATALOG #



Screwdriver Bit Holder 4100-023
NAME CATALOG #

Screwdriver Bit- Single Slot
 Trinkle Shank 4100-024
NAME CATALOG #

Screwdriver Bit- Cruciate
 Trinkle Shank 4100-026
NAME CATALOG #

Screwdriver Bit-Phillips
 Trinkle Shank 4100-025
NAME CATALOG #

Screwdriver Bit-A.S.I.F.
 Trinkle Shank 4100-027
NAME CATALOG #

NOTE: The 4100-023 bit holder fits only the 4100-027 A.S.I.F. screwdriver bit.

PULSE LAVAGE



Pulse Lavage Coupler 7100-474
NAME CATALOG #
 NOTE: Uses 4740 tubing sets

