







Continuous Carbon Fibers Reinforced Polymer

- Optimal Modulus of Elasticity
- Unparalleled Fatigue Strength
- Radiolucent
- MRI/CT- Artifact Free
- Essential in Radiation Therapy
- Plates: No Cold - Welding
- Common Sizes & Surgical Technique

Piccolo Composite® Nails & Plates




Clearly Better

www.carbo-fix.com




The Piccolo Composite® Portfolio





Unique Biomechanical & Imaging Properties

*Clear visibility
through the plate*




Distal Fibula Plate

*Clear visibility
through the nail*




Trochanteric (PF) Nail

*MRI & CT: No artifacts
Clear image of adjacent tissue*




Humeral Nail

Improved modulus of elasticity




Alok D Sharan, Simon Y Tang, Alexander R. Basic Science of Spinal Diseases
Japsee brothers medical publishers.2013; 159-160

Superior bending strength



* Performed on 4.5mm Diaphyseal Plates

Unparalleled fatigue strength



* Bending test: 1800N at 5Hz on Ø 8.5mm rods

Plates

Implant	Number of Holes (Length, mm)	Orientation	Thickness (mm)
Distal Radius-Standard	3, 4, 7 (52, 60, 90)	Right or Left	2.4
Distal Radius-Narrow	3, 4 (52, 60)	Right or Left	2.4
Distal Radius-Triangular	3, 4 (54, 63)	Right or Left	2.4
Proximal Humerus	3, 4*, 6, 8, 12* (102, 120*, 156, 192, 262*)	Bi-Lateral	3.7
Distal Fibula	3*, 4, 5*, 6, 9, 11 (72*, 85, 96*, 114, 155, 184)	Right or Left	2.3
One Third Tubular	5, 6, 7, 9, 11 (56, 67, 78, 98, 122)	Bi-Lateral	2.3
Diaphyseal-Narrow	7, 9 (154, 190)	Bi-Lateral	4.7
Diaphyseal-Broad	9, 11, 13 (160, 185, 220)	Bi-Lateral	5.3
Diaphyseal No Holes	(160, 185, 220)	Bi-Lateral	5.3
Distal Femur	6, 8, 10, 12, 14, 16 (171, 207, 243, 279, 315, 350)	Right or Left	5.1
MTP	Dorsiflexion: 0°, 5° (50)	Right or Left	2.0

Nails

Implant	Length (mm)	Diameter (mm)	Prox. Diameter (mm)
Humerus	180, 200, 220, 240, 280	8.5	10.0
Proximal Humerus	150	8.0	11.0
Tibia	260, 280, 300, 320, 340, 360, 380, 400	10.0, 11.0	11.5
Femur	300, 320, 340, 360, 380, 400, 420	10.0, 11.0, 12.0	11.5, 11.5, 12.0
Ankle Arthrodesis	160, 200, 240	10.0, 12.0	12.0

Trochanteric (PF) Nail	180, 300, 320, 340, 360, 380, 400, 420, 440*, 460*	11.0	17.0
Lag Screw	80, 85, 90, 95, 100, 105, 110	10.4	---

* Special order

MANUFACTURED BY:

CarboFix Orthopedics Ltd.
11 Ha'hoshlim St.,
Herzeliya 4672411, Israel
Tel: +972-9-9511511
Fax: +972-9-9548939
E-Mail: info@carbo-fix.com

U.S.A. OFFICE:

CarboFix Orthopedics Inc.
7183 Beach Dr SW Ste 1
Ocean Isle Beach, NC 28469, USA
Tel: 1-800-408-0120
Fax: 877-705 3567
E-Mail: usa@carbo-fix.com

EC AUTHORIZED REPRESENTATIVE:

MEDNET GmbH
Borkstrasse 10, 48163 Münster
Germany



OUDSHOORN
chirurgische techniek

A Wethouder van Eschstraat 38, 5342 AT OSS
T 0412 - 63 24 00
W www.oudshoornbv.com
E info@oudshoornbv.com

Patents are pending

Caution: In the U.S.A., federal law restricts this device to sale by or on the order of a physician.