


INTEOS[®] HCS(e) 3.5 Clavicle Plate sup. (ext.)
Anatomically optimized clavicle plate system



Advantages of the INTEOS® HCS(e) 3.5 Clavicle Plate sup. (ext.)

Die INTEOS® HCS(e) 3.5 Clavicle Plate sup. (ext.) enable a stable treatment of diaphyseal clavicle fractures due to their anatomically correct fit and on the proven basis of the **angle-stable & polyaxial** plate-screw connection from Hofer Medical Solutions.

- + Precise anatomical pre-bending of the plate
- + Available in two radii of curvature, with standard or extended geometry
- + Up to 3 different plate lengths in left and right versions for optimal adaptation to the patient's anatomy
- + Fully retractable screw heads with high polyaxial angular stability
- + Minimized soft tissue irritation thanks to the rounded geometry of the plate and the deep-seated screw heads
- + Angle-stable and polyaxial locking mechanism
- + INTEOS® 3.5 compression holes for standard screws in the middle section of the plates to ensure highest stability
- + System added value: compatible with screws and instruments of the INTEOS® 3.5 System



L

INTEOS® HCS 3.5 Clavicle Plate sup. as av



Plate thickness 3.0 mm
Compatible to:
Shaft area: INTEOS® 3.5 screws / ForceDRIVE 10

Titanium	Number of Holes	Length	Side	PU
750-100-060-007	5/2	92 mm	left	1
750-100-060-009	7/2	110 mm	left	1
750-100-060-011	9/2	127 mm	left	1
750-100-061-007	5/2	92 mm	right	1
750-100-061-009	7/2	110 mm	right	1
750-100-061-011	9/2	127 mm	right	1

Additional options

S	xxx-xxx-xxx-xxx-S	Sterile	
T	xxx-xxx-xxx-xxx-T	Trial	Do not implant!

L

INTEOS® HCSe 3.5 Clavicle Plate ext. sup. as av




Plate thickness 3.0 mm
Compatible to:
Shaft area: INTEOS® 3.5 screws / ForceDRIVE 10

Titanium	Number of Holes	Length	Side	PU
750-100-060-107	5/2	89 mm	left	1
750-100-060-109	7/2	107 mm	left	1
750-100-060-111	9/2	125 mm	left	1
750-100-061-107	5/2	89 mm	right	1
750-100-061-109	7/2	107 mm	right	1
750-100-061-111	9/2	125 mm	right	1

Additional options

S	xxx-xxx-xxx-xxx-S	Sterile	
T	xxx-xxx-xxx-xxx-T	Trial	Do not implant!





AUSTRIA

HOFER GMBH & CO KG
Jahnstrasse 10-12
A 8280 Fürstenfeld
TEL: +43(0)3382/53388
FAX: +43(0)3382/53093
office@hofer-medical.com
www.hofer-medical.com
FN21826y
UID ATU30764704



GERMANY

HOFER-MEDICAL GMBH
Ruhrtalstraße 52-60
D 45239 Essen
TEL: 0800-700 80 22 (innerhalb BRD)
FAX: +49 201 56599 1069
office@hofer-medical.de
www.hofer-medical.com
HRB24645
UID DE290596924



ITALY

HOFER MEDICAL ITALIA S.R.L.
Via della Rena, 26
I - 39100 Bolzano

Tel: +39 393 8950481
FAX: +39 06 81151148
office@hofer-medical.it
www.hofer-medical.com
RM-1454597
UID IT13534721009



POLAND

HOFER GMBH & CO KG SPÓŁKA
KOMANDYTOWA ODDZIAŁ W POLSCE
Radkowska 9
PL 57-402 Nowa Ruda

Tel.Kom.:+48603224112
Fax: +48 748712273
office@hofer-medical.pl
www.hofer-medical.com
REGON 141527293
UID PL1070011771

creating partnership

For us, a partnership is the result of a long and happy relationship in all areas of our work. Reliable bone healing for patients, optimal handling of the systems for surgeons and surgical personnel, as well as the simple preparation of the instruments, constantly strengthen this partnership.

For us, "creating" means more than just finding solutions in the form of innovative products.

It is the result of high standards, constant development, innovative products and excellent service for patients, surgeons and surgical personnel around the clock.

Please do not hesitate to contact us if you have questions about our company, our employees or our production methods.

Your Hofer Medical SolutionsTeam

vCard

