

INTEOS[®] 2.5 Coronoid Plate

Innovative osteosynthesis system for coronoid fractures



Advantages of the INTEOS® 2.5 Coronoid Plate

The new INTEOS® 2.5 coronoid plates from Hofer Medical Solutions are characterized by their anatomically adapted shape and the extra flat design while maintaining the maximum possible stability.

- + Extra flat and rounded design to minimize irritation
- + Fully retractable screw heads with high polyaxial angular stability
- + Bendable flap for support and fixation of the sublimus tubercle to restore elbow stability
- + Subchondral screw placement is possible
- + The coronoid plates are anatomically pre-bent, but can be bent in exceptional cases
- + K-wire holes for temporary plate and fragment fixation
- + A hook is located proximally to fix small fragments outside the screw range
- + System added value: compatible with screws and instruments of the INTEOS® 2.5 System



S

INTEOS® 2.5 Coronoid Plate as av T5/1



Plate thickness 1,5 mm
Compatible to INTEOS® 2.5 screws / ForceDRIVE 6

Titanium	Number of holes	Side	Length	PU
777-170-000-006	3 Head / 3 Shaft	left	38 mm	1
777-170-001-006	3 Head / 3 Shaft	right	38 mm	1



Additional options

S xxx-xxx-xxx-xxx-**S** Sterile





AUSTRIA

HOFER GMBH & CO KG
Jahnstrasse 10-12
A 8280 Fürstenfeld
TEL: +43(0)3382/53388
FAX: +43(0)3382/53093
office@hofer-medical.com
www.hofer-medical.com
FN21826y
UID ATU30764704



GERMANY

HOFER-MEDICAL GMBH
Ruhrtalstraße 52-60
D 45239 Essen
TEL: 0800-700 80 22 (innerhalb BRD)
FAX: +49 201 56599 1069
office@hofer-medical.de
www.hofer-medical.com
HRB24645
UID DE290596924



ITALY

HOFER MEDICAL ITALIA S.R.L.
Via della Rena, 26
I - 39100 Bolzano

Tel: +39 393 8950481
FAX: +39 06 81151148
office@hofer-medical.it
www.hofer-medical.com
RM-1454597
UID IT13534721009



POLAND

HOFER GMBH & CO KG SPÓŁKA
KOMANDYTOWA ODDZIAŁ W POLSCE
Radkowska 9
PL 57-402 Nowa Ruda

Tel.Kom.:+48603224112
Fax: +48 748712273
office@hofer-medical.pl
www.hofer-medical.com
REGON 141527293
UID PL1070011771

creating partnership

For us, a partnership is the result of a long and happy relationship in all areas of our work . Reliable bone healing for patients, optimal handling of the systems for surgeons and surgical personnel, as well as the simple preparation of the instruments, constantly strengthen this partnership.

For us, "creating" means more than just finding solutions in the form of innovative products.

It is the result of high standards, constant development, innovative products and excellent service for patients, surgeons and surgical personnel around the clock.

Please do not hesitate to contact us if you have questions about our company, our employees or our production methods.

Your Hofer Medical SolutionsTeam

vCard

