



bee-system II
Li-Ion - Technology



sma bo
small bone

Contact us for further technical information or
Make an appointment for a live demonstration.

sales@medical-bees.de


We appreciate your interest.

MB07-10-001-00E © Changes in technical specifications and dimensions are at our discretion. Subject to modification. Without our explicit written approval this print media may not be copied or reprinted as a whole or in extracts.

medicalbees
System-driven Surgery
Operieren mit System

Powertools

sma bo
small bone



The innovation for small bone surgery

The Number 1 in



Weight: 33% lighter than beesystemII



Design: 1 universal handpiece, various attachments with **One-Click adaptation** for drill, saw and pin-driver



Technology: Smabo with smart charging technology via power cord supply and easy magnetic connection

The proven concept of a non-sterilizable Powerpack is the heart of the system for a unique durability. All integrated components: motor, battery and electronics in removable unit.



F-20-100-00	smabo universal handpiece	F-20-200-50	smabo Drill attachment AO-small
F-20-900-00	smabo Powerpack	F-20-200-60	smabo Jacobs attachment with key
F-20-900-20	smabo sterile funnel	F-20-300-00	smabo oscillating saw attachment with key
F-20-900-40	smabo charging cable with power supply for smabo Powerpack F-20-900-00	F-20-500-00	smabo K-wire/ pin-driver
F-20-000-03	smabo presentation case	F-20-910-00	smabo tray with lid